

# birth pool inflation guide

- Here are the basic steps to take when inflating your birth pool.
- Firstly, remove the pool from packing and warm to room temperature for 72 hours before inflation
- Using the electric air pump provided, select the smallest valve adapter
- Inflate the floor first
- Unscrew cap on side valve (top one). Select the appropriate valve adapter on the electric air pump, start inflating at the bottom ring
- As soon as the pool ring is firm, stop the electric pump and replace side valve cover
- Repeat this process for both bottom and middle valve
- When inflating the top ring, leave it slightly deflated to allow the liner to be put on
- Inflate the seat
- Put on pool liner
- After the liner is in place, inflate the top ring fully and when the pool is firm, screw on the side valve cap.
- Fill the pool with warm water
- That's it you're all ready to go!

